


Product Name: Recombinant Human DZIP3 Protein
Catalog #: PEH8755



Summary

Name	DZIP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86Y13
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPVPEIIQRM LSCYQQGIAL QSITGSQRIE IEELQNEEEE LSPPLMEYNI NVKSHPEIQF AEINKDGTSI PSESSTESLK DLQEVKSKQR KKKKTKNKKK KDSKEDQVPY VVEKEEQLRK EQANPHSVSR LIKDDASDVQ EDSAMEDKFY SLDELHILDM

image



Background

Alternative Names	E3 ubiquitin-protein ligase DZIP3 (EC 2.3.2.27) (DAZ-interacting protein 3) (RING-type E3 ubiquitin transferase DZIP3) (RNA-binding ubiquitin ligase of 138 kDa) (hRUL138)
--------------------------	--

Note

Product Name: Recombinant Human DZIP3 Protein
Catalog #: PEH8755



For Research Use Only , Not for Diagnostic Use