


Summary

Name	EIPR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53HC9
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IISLADNHIL LWDLQESSQ AVLASSASLE GKGQLKFTSG RWSPHHNCTQ VATANDTTLR GWDTRSMSQI YCIENAHGQL VRDLDFNPNK QYYLASCDD CKVKFWDTRN VTEPVKLEE HSHWVWNVRY NHSHDQLVLT GSSDSRVILS NMVSISSPEF GHLVDDDDIS DQEDHRSEEK SKEPLQDNVI

image



Background

Alternative Names	EARP and GARP complex-interacting protein 1 (Endosome-associated recycling protein-interacting protein) (Golgi-associated retrograde protein-interacting protein) (Tumor-suppressing STF cDNA 1 protein) (Tumor-suppressing subchromosomal transferable fragment candidate gene 1 protein)
--------------------------	--

Product Name: Recombinant Human EIPR1 Protein
Catalog #: PEH8756



Note

For Research Use Only , Not for Diagnostic Use