


Product Name: Recombinant Human METTL13 Protein
Catalog #: PEH8757



Summary

Name	EFNMT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N6R0
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVPFLSVGGD IGVRTVQHGD CSPLSGDYVI EDVQGDDKRY FRRLIFLSNR NVVQSEARLL KDVSHKAQKK RKKDRKKQRP ADAEDLPAAP

image



Background

Alternative Names	eEF1A lysine and N-terminal methyltransferase (eEF1A-KNMT) (Methyltransferase-like protein 13) [Includes: eEF1A lysine methyltransferase (EC 2.1.1.-); eEF1A N-terminal methyltransferase (EC 2.1.1.-)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human METTL13 Protein
Catalog #: PEH8757

