


Product Name: Recombinant Human EFHD1 Protein
Catalog #: PEH8759



Summary

Name	EFHD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUP0
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERRLRREEAE ESGPQLAPLG APAPEPKPEP EPPARAPTAS ADAELSAQLS RRLDINEGAA RPRRCRVFNP YTEFPEFSRR LIKDLESMFK LYDAGRDFGI DLMELKLMME KLGAPQTHLG LKSMIKEVDE DFDGKLSFRE FLLIFHKAAA GELQEDSGLM ALAKLSEIDV ALEGVKGAKN FFEAKVQALS SASKFEAEK AEQDERKREE

image



Background

Alternative Names	EF-hand domain-containing protein D1 (EF-hand domain-containing protein 1) (Swiprosin-2)
--------------------------	--

Note

Product Name: Recombinant Human EFHD1 Protein
Catalog #: PEH8759



For Research Use Only , Not for Diagnostic Use