


Product Name: Recombinant Human ELMO3 Protein
Catalog #: PEH8760



Summary

Name	ELMO3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BJ8
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLRLRQTERL HQEGLAPPI LELREKLKPE LMGLIRQQRL LRLCEGTLFR KISSRRRQDK LWFCCLSPNH KLLQYGDMEE GASPPTLESL PEQLPVADMR ALLTGKDCPH VREKGSQKQN KDLYELAFSI SYDRGEEAY LNFIAPSKRE FYLWTDGLSA LLGSPMGSEQ TRLDLEQLLT METKLRLLEL ENVPIPERPP PVPPPPTNFN FCYDCSIAEP

image



Background

Alternative Names Engulfment and cell motility protein 3

Note

Product Name: Recombinant Human ELMO3 Protein
Catalog #: PEH8760



For Research Use Only , Not for Diagnostic Use