


Product Name: Recombinant Human EVI5L Protein
Catalog #: PEH8764



Summary

Name	EVI5L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96CN4
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKRLRTENRL LKQRIETLEK GQVTRAQAE ENYVIKRELA VVRQQCSSAA EDLQKAQSTI RQLQEQQENP RLTEDFVSHL ETELEQSRLR ETETLGALRE MQDKVLDMEK RNSSLPDENN VAQLQEELKA LKVREGQAVA STRELKLQLQ ELSDTWQAHL ARGGRWKESP RKLVGELQD ELMSVRLREA QALAEGRRLR QRVVELETQD HIHRNLLNRV EAERAALQEK LQYLAAQNKG LQTQLSESRR KQAEAECKSK EEVMAVRLRE ADSMAAVAEM RQRIAELEIQ REEGRIQGQL NHSDSSQYIR ELKDQIEELK

image



Background

Alternative Names EVI5-like protein (Ecotropic viral integration site 5-like protein)

Product Name: Recombinant Human EVI5L Protein
Catalog #: PEH8764



Note

For Research Use Only , Not for Diagnostic Use