



Product Name: Recombinant Human EXOC3L4 Protein
Catalog #: PEH8765



Summary

Name	EX3L4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q17RC7
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PELQSPKEAE EPQTPAQGSR RTSSRKEPNA HRKDGTRLGL GSLRQAFSRA SQRALTQVSK EDTGLFRRSS CSLFRSFRQA LNDGPATGHS QATPEVPSGV MNGVSQQAST GAASEELKPE AEGKSVADLI TERQLLAFFE QLLRLETLV AEKASRTFEQ DPTAFARRAM DVCLLYDGLA AEIGAIVRET LDSDGVDAAA LAELARVSA EEEAHPSPPD DGDFLRTPRR WRQHWEEAVR RSAQERVRRP GAGWAFGEAE GASGLAQLLA ELGGLVRRDL QKVRQEVQPA YAAAGFPAWE VYLRAFHSVA AQRLQELARD ARGCEQLYL

image



Background

Alternative Names Exocyst complex component 3-like protein 4

Product Name: Recombinant Human EXOC3L4 Protein
Catalog #: PEH8765



Note

For Research Use Only , Not for Diagnostic Use