


Product Name: Recombinant Human FCHSD2 Protein
Catalog #: PEH8766



Summary

Name	FCSD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94868
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVFHKQPFQ FQPCSDTSR QLESETGTTE EHSLNKEARK WATRVAREHK NIVHQQRVLN DLECHGAAVS EQSRAELEQK IDEARENIRK AEIILKAEA RLDLLKQIGV SVDTWLKSAM NQVMEELENE RWARPPAVTS NGTLHSLNAD TEREEGEEFE DNMDVFDDSS SSPSGTLRNY PLTCKVVVYSY KASQPDELTI EEHEVLEVIE DGD MEDWVKA RNKVGQVGYV PEKYLQFPTS

image



Background

Alternative Names F-BAR and double SH3 domains protein 2 (Carom) (Protein nervous wreck 1) (NWK1) (SH3 multiple domains protein 3)

Note

Product Name: Recombinant Human FCHSD2 Protein
Catalog #: PEH8766



For Research Use Only , Not for Diagnostic Use