


Product Name: Recombinant Human FCHO2 Protein
Catalog #: PEH8767



Summary

Name	FCHO2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0JRZ9
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WGEKNSGFDV LYHNMKHGQI STKELADFVR ERATIEEAYS RSMTKLAKSA SNYSQLGTFA PVWDVFKTST EKLANCHLDL VRKLQELIKE VQKYGEEQVK SHKKTKEEVA GTLEAVQTIQ SITQALQKSK ENYNAKCVEQ ERLKKEGATQ REIEKAAVKS KKATDITYKLY VEKYALAKAD FEQKMTETAQ KFDIEETHL IHIKEIIGSL SNAIKEIHLQ IGQVHEEFIN NMANTTVESL IQKFAESKGT GKERPGLIEF

image



Background

Alternative Names F-BAR domain only protein 2

Product Name: Recombinant Human FCHO2 Protein
Catalog #: PEH8767



Note

For Research Use Only , Not for Diagnostic Use