


Product Name: Recombinant Human FBXL7 Protein
Catalog #: PEH8771



Summary

Name	FBXL7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJT9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEGKGSISIS SDVSSSTDHT PTKAQKNVAT SEDSDLSMRT LSTPSPALIC PPNLPGFQNG RGSSTSSSSI TGETVAMVHS

image



Background

Alternative Names F-box/LRR-repeat protein 7 (F-box and leucine-rich repeat protein 7) (F-box protein FBL6/FBL7)

Note

For Research Use Only , Not for Diagnostic Use