


Product Name: Recombinant Human FRMD4A Protein
Catalog #: PEH8772



Summary

Name	FRM4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2Q2
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSAKKDMLAA LKSRQEAL EE TLRQRLEELK KLCLREAELT GKLPVEYPLD PGEEPPIVRR RIGTAFKLDE QKILPKGEEA ELERLEREFA IQSQITEAAR RLASDPNVSK KLKKQRKTSY LNALKKLQEI ENAINENRIK SGKKPTQRAS LIIDDGNIA S EDSSLSDALV LEDEDSQVTS TISPLHSPHK GLPPRPPSHN RPPPPQSLEG LRQMHYHRND YDKSPIKPKM WSESSLDEPY

image



Background

Alternative Names FERM domain-containing protein 4A

Note

Product Name: Recombinant Human FRMD4A Protein
Catalog #: PEH8772



For Research Use Only , Not for Diagnostic Use