


Product Name: Recombinant Human FPGT Protein
Catalog #: PEH8773



Summary

Name	FPGT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14772
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQELAYNQQL SEKLKRKELP LGVQYHVFVD PAGAKIGNGG STLALQCLE KLYGDKWNSF TILLHSGGY SQRLPNASAL GKIFTALPLG NPIYQMLELK LAMYIDFPLN MNPGILVTCA DDIELYSIGE FEFIRFDKPG FTALHPSSL TIGTTHGVFV

image



Background

Alternative Names	Fucose-1-phosphate guanylyltransferase (EC 2.7.7.30) (GDP-L-fucose diphosphorylase) (GDP-L-fucose pyrophosphorylase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FPGT Protein
Catalog #: PEH8773

