



Product Name: Recombinant Human FGD6 Protein
Catalog #: PEH8774



Summary

Name	FGD6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZV73
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSMEDADANV YEVEEPEAP DGQLQLGPRH QHSSSGASQE EQNDLGLGDL PSDEEEIINS SDEDDVSES SKGEPDPLED KQDEDNGMKS KVHHIAKEIM SSEKVFVDVL KLLHIDFRDA VAHASRQLGK PVIEDRILNQ ILYYLPQLYE LNRDLLKELE ERMLHWTEQQ RIADIFVKKG PYLKMYSTYI KEFDKNIALL DEQCKKNPGF AAVVREFEMS PRCANLALKH YLLKPVQRIP QYRLLLTGYL KNLIEDAGDY RDTQDALAVV IEVANHANDT

image



Background

Alternative Names	FYVE, RhoGEF and PH domain-containing protein 6 (Zinc finger FYVE domain-containing protein 24)
--------------------------	---

Product Name: Recombinant Human FGD6 Protein
Catalog #: PEH8774



Note

For Research Use Only , Not for Diagnostic Use