


Product Name: Recombinant Human TUBGCP4 Protein
Catalog #: PEH8777



Summary

Name	GCP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UGJ1
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEFFIKQGPS SGNVSAQPEE DEEDLGIGGL TGKQLRELQD LRLIEENML APSLKQFSLR VEILPSYIPV RVAEKILFVG ESVQMFENQN VNLTRKGSIL KNQEDTFAAE LHRLKQQPLF SLVDFEQVVD

image



Background

Alternative Names	Gamma-tubulin complex component 4 (GCP-4) (hGCP4) (Gamma-ring complex protein 76 kDa) (h76p) (hGrip76)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use