


**Product Name: Recombinant Human GATB Protein**  
**Catalog #: PEH8789**



## Summary

<b>Name</b>	GATB_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O75879
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	36 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	NRAGVGLLEV VLEPDMSCGE EAATAVRELQ LILQALGTSQ ANMAEGQLRV DANISVHHPG EPLGVRTEVK NLNSIRFLAK AIDYEIQRQI NELENGGEIL NETRSFHHKL GCTMSMRDKE GKQDYRFMPE PNLPLVLYD ATSLPAGADP QQVINIDQIR ETLPELPSVT REKLVQQYGM LLEHSFTLLN EVGLLEFFQN VIKETRAEPK KVTSWVLNTF LGYLKQQNLA VSESPVTPSA LAELDLLDS RTISSAAKQ VFEELWKREG KTPGQIVSEK QLELMQDQGA LEQLCHSVME

## image



## Background

<b>Alternative Names</b>	Glutamyl-tRNA(Gln) amidotransferase subunit B, mitochondrial (Glu-AdT subunit B) (EC 6.3.5.-) (Cytochrome c oxidase assembly factor PET112 homolog)
--------------------------	---

**Product Name: Recombinant Human GATB Protein**  
**Catalog #: PEH8789**

---



**Note**

For Research Use Only , Not for Diagnostic Use