

Product Name: Recombinant Human GOLGA8CP Protein
Catalog #: PEH8794



Summary

Name	GOG8C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NN73
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQMAEPLPPE PPAVPSEVEL QHLRKELERV AGALQAQVEY NQRISLLNEG QKERLREQEE RLQEQQERLR EQEERLQQLA EPQNSFKELN NENKSVLQLE QQVKELQEKL GKRLAHPVAS AQKEPEAAVP APGPGGESSG FMDHLEEKAD LSELVEKEEL GFFQYYRERC HQKVYHPITK PGGSAKDAAP GGGHHQAGPG QGGDEGEAAG AAGDGVAAGG DYKGHSKFLV

image



Background

Alternative Names Golgin subfamily A member 8C

Note

Product Name: Recombinant Human GOLGA8CP Protein
Catalog #: PEH8794



For Research Use Only , Not for Diagnostic Use