

Product Name: Recombinant Human HAUS7 Protein
Catalog #: PEH8799



Summary

Name	HAUS7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99871
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QWASASAKSE EEEKLAELAR QLQESAACLH ALRTEYFAQH EQGAAAGAAD ISTLDQKLRL VTSDFHQLIL AFLQVYDDEL GECCQRPGPD

image



Background

Alternative Names HAUS augmin-like complex subunit 7 (26S proteasome-associated UCH37-interacting protein 1) (UCLH5-interacting protein) (X-linked protein STS1769)

Note

For Research Use Only , Not for Diagnostic Use