


Product Name: Recombinant Human HEATR3 Protein
Catalog #: PEH8800



Summary

Name	HEAT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z4Q2
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPQFSPTGDC QAEAAAAANG TGGEEDDGPA AELLEKLQHP SAEVRECACA GLARLVQQRP ALPGLARRDA VRRLGPLLLD PSLAVRETAA GALRNLSACG GFEVCDDMVT KDIMTPLVAL LKECSAGLDS NEMSLQEKKD QNRNSIENIA NETVNVLWNI CECSSRAVSI FNKEGCLEIV LKYLRFPTN VDLAISVAYC LQTVTEDNPE LLKSFSATAL NMLESALLSP VSSMESLLK TLVAGTIWNL KDIIPCKSQA EIINALLKIL SEVLGMDAGE MVIQMKEAET QRLKTAEEAE EILENTNGDD LIEDEMEGI SHKRRVRRKT FVSDLLPPTD

image



Background

Alternative Names	HEAT repeat-containing protein 3
--------------------------	----------------------------------

Product Name: Recombinant Human HEATR3 Protein
Catalog #: PEH8800



Note

For Research Use Only , Not for Diagnostic Use