


**Product Name: Recombinant Human HNRNPLL Protein**  
**Catalog #: PEH8802**



## Summary

|                              |                                                                                                                                                                                                                                                                                            |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | HNRLH_HUMAN                                                                                                                                                                                                                                                                                |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                        |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                           |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                      |
| <b>Accession #</b>           | Q8WVV9                                                                                                                                                                                                                                                                                     |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                     |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                      |
| <b>Molecular Mass</b>        | 30 kDa                                                                                                                                                                                                                                                                                     |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                         |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                             |
| <b>Sequence</b>              | TYEEDREYES QAKRLKTEEG EIDYSAEEGE NRREATPRGG GDGGGGGRSF<br>SQPEAGGSHH KVSVPVHV RGLCESVVEA DLVEALEKFG TICVMMMPF<br>KRQALVEFEN IDSAKECVTF AADEPVYIAG QQAFFNYSTS KRITRPGNTD<br>DPSGGNKVLL LSIQNPLYPI TVDVLYTCN PVGKVQRIVI FKRNGIQAMV<br>EFESVLCAQK AKAALNGADI YAGCCTLKIE YARPTRLNVI RNDNDSWDYT |

## image



## Background

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Heterogeneous nuclear ribonucleoprotein L-like (hnRNPLL) (Stromal RNA-regulating factor) |
|--------------------------|------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human HNRNPLL Protein**  
**Catalog #: PEH8802**



---

For Research Use Only , Not for Diagnostic Use