


Product Name: Recombinant Human GGT2P protein
Catalog #: PEH8806



Summary

Name	GGT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P36268
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRLADTYEML AIEGAQAFYN GSLMAQIVKD IQAAGGIVTA EDLNNYRAEL IEHPLNISLG DAVLYMPSAR LSGPVLALIL NILKGYNFSR ESVETPEQKG LTYHRIVEAF RFAYAKRTLL GDPKFVDVTE VVRNMTSEFF AAQLRSQISD HTTHPISYYK PEFYTPDDGG TAHLSVVAED GSAVSATSTI NLYFGSKVCS PVSGILFNNE WTTSALPAFT NEFGAPPSPA NFIQPGKQPL LSMCLTIMVG QDGQVRMVG AAGGTQITTD TALAIYNLW FGYDVKRAVE EPRLHNKLLP NVTTVERNID

image



Background

Alternative Names Inactive glutathione hydrolase 2 (Gamma-glutamyltransferase 2) (Inactive gamma-glutamyltranspeptidase 2) (GGT 2)

Product Name: Recombinant Human GGT2P protein
Catalog #: PEH8806



Note

For Research Use Only , Not for Diagnostic Use