



## Summary

<b>Name</b>	UBP54_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q70EL1
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	43 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TVIYNVENDS MSQSSRDTGH LTDSECNQKH TSKKGLIER KRSSGRVRRK GDEPQASGYH SEGETLKEKQ APRNASKPSS STNRLRDFKE TVSNMIHNRP SLASQTNVGS HCRGRGGDQP DKKPPRTLPL HSRDWEIEST SSESKSSSSS KYRPTWRPKR ESLNIDSIFS KDKRKHCGYT QLSPFSEDSA KEFIPDEPSK PPSYDIKFGG PSPQYKRWGP ARPGSHLLEQ HPRLIQRMES GYESSERNSS SPVSLDAALP ESSNVYRDPS AKRSAGLVPS WRHIPKSHSS SILEVDSTAS MGGWTKSQPF SGEEISSKSE LDELQEEVAR RAQEQLRRK REKELEAAKG FNPHPSRFMD

## image



## Background

**Alternative Names** Inactive ubiquitin carboxyl-terminal hydrolase 54 (Inactive ubiquitin-specific

**Product Name: Recombinant Human USP54 Protein**  
**Catalog #: PEH8808**



---

peptidase 54)

**Note**

For Research Use Only , Not for Diagnostic Use