


Product Name: Recombinant Human IFT46 Protein
Catalog #: PEH8811



Summary

Name	IFT46_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQC8
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EENNKEKKKT SQLTPQRGFS ENEDDDDDDD DSSETDSDS DDDEEHGAPL EGAYDPADYE HLPVSAEIKE LFQYISRYTP QLIDLHKLK PFIPDFIPAV GDIDAFKVP RPDGKPDNLG LLVLDEPSTK QSDPTVLSLW LTENSKQHNI TQHMKVKSLE DAEKNPKAID TWIESISELH RSKPPATVHY TRPMPDIDTL MQEWSPEFEE LLGKVSPLTA EIDCSLAEYI DMICAILDIP VYKSRIQSLH LLFSLYSEFK NSQHFALAE

image



Background

Alternative Names Intraflagellar transport protein 46 homolog

Product Name: Recombinant Human IFT46 Protein
Catalog #: PEH8811



Note

For Research Use Only , Not for Diagnostic Use