


Product Name: Recombinant Human IQCE Protein
Catalog #: PEH8813



Summary

Name	IQCE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IPM2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTGAKRQKKM GSALLSLRS VQELTEENQS LKEDLDRVLS TSPTISKTQG YVEWSKPRLL RRIVELEKKL SVMESSKSHA AEPVRSHPPA CLASSSALHR QPRGDRNKDH ERLRGAVRDL KEERTALQEQ LLQRDLEVKQ LLQAKADLEK ELECAREGEE ERREREEVLR EEIQTLSKL QELQEMKKEE KEDCPEVPHK AQELPAPTPS SRHCEQDWPP

image



Background

Alternative Names IQ domain-containing protein E

Note

Product Name: Recombinant Human IQCE Protein
Catalog #: PEH8813



For Research Use Only , Not for Diagnostic Use