


Product Name: Recombinant Human ACAD8 Protein
Catalog #: PEH8814



Summary

| | |
|------------------------------|---|
| Name | ACAD8_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9UKU7 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 26 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | LNEEQKEFQK VAFDFAAREM APNMAEWDQK ELFPVDVMRK AAQLGFGGVY IQTDVGGGSL SRLDTSVIFE ALATGCTSTT AYISIHNMCA WMIDSFGNEE QRHKFCPPLC TMEKFASYCL TEPGSGSDAA SLLTSAKKQG DHYILNGSKA FISGAGESDI YVVMCRTGGP GPKGISCIVV EKGTPGLSFG KKEKKVWNS QPTRAVIFED CAVPVANRIG SEGQGFLIAV |

image



Background

| | |
|--------------------------|---|
| Alternative Names | Isobutyryl-CoA dehydrogenase, mitochondrial (IBDH) (EC 1.3.8.-) (Activator-recruited cofactor 42 kDa component) (ARC42) (Acyl-CoA dehydrogenase family member 8) (ACAD-8) |
|--------------------------|---|

Product Name: Recombinant Human ACAD8 Protein
Catalog #: PEH8814



Note

For Research Use Only , Not for Diagnostic Use