


Product Name: Recombinant Human JHY Protein
Catalog #: PEH8815



Summary

Name	JHY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NUN7
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAQQLKNYQE HWSQYESTKS SNVPRGQPSD MVNDHQPSRR PAKLKIRKQC KHQNLKSST TEEVTASQGN QNNPPRQQQN QNKPLDTSTK PESIVIMHAS NNDVQASRAL RSHNLKETS N TFAPPKQAFD KVLSTKSTGC DSGLNVNKER GHKDQEEKRF SYQQLHTLSD MDLNNLNELS KRHVLLSQKG

image



Background

Alternative Names Jhy protein homolog

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human JHY Protein
Catalog #: PEH8815

