


Product Name: Recombinant Human KLHL25 Protein
Catalog #: PEH8819



Summary

Name	KLH25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0H3
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RTKLIMDEAL RCKTRILQND GVVTSPCARP RKAGHTLLIL GGQTFMCDKI YQVDHKAKEI IPKADLPSPR KEFSASAIGC KVVVTGGRGS ENGVSKDVVV YDTVHEEWSK AAPMLIARFG HGSAELENCL YVVGHTSLA GVFPASPSVS LKQVEKYDPG ANKWMMVAPL RDGVSNAAVV SAKLKLFFVG GTSIHRDMVS KVQCYDPSEN RWTIKAACPQ PWRYTAAAVL GSQIFIMGGD

image



Background

Alternative Names Kelch-like protein 25 (Ectoderm-neural cortex protein 2) (ENC-2)

Note

Product Name: Recombinant Human KLHL25 Protein
Catalog #: PEH8819



For Research Use Only , Not for Diagnostic Use