


Product Name: Recombinant Human KNDC1 Protein
Catalog #: PEH8821



Summary

Name	KNDC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q76NI1
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFCDLYWDEK LLQNLFKVVN GQASPPSTA EEAGSQLEGS QSPRSPSSKR PSLHRLGKEK PAMARTSSRA PCSPTSVDV DSDALSRGNF EVGFRPQRSV KAERAQQPEA GEDRRPAGGA SDVEAVTRLA

image



Background

Alternative Names Kinase non-catalytic C-lobe domain-containing protein 1 (KIND domain-containing protein 1) (Cerebral protein 9) (Protein very KIND) (v-KIND) (Ras-GEF domain-containing family member 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KNDC1 Protein
Catalog #: PEH8821

