

Product Name: Recombinant Human KRBA2 Protein
Catalog #: PEH8826



Summary

Name	KRBA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZNG9
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WQQPMLFQEA TAFENMTKDW NYLEGSQKDC YRDTMLDSYE NTPVQGSFLQ LSMMPQRAGN DPPGVSNASE MEMEISNMRE KFLMSVTKLV ESKSYNSKVF SKEYFQTIK EVKEAKEKGK KSSRDYRRAA KYDVISVQGT EKLIEATHGE RDRIYYVHK EELFDILHDT HLSIGHGGRT RMLKELQGKY GNVTKEVIVL YLTLCQCHQ

image



Background

Alternative Names KRAB-A domain-containing protein 2

Note

Product Name: Recombinant Human KRBA2 Protein
Catalog #: PEH8826



For Research Use Only , Not for Diagnostic Use