


Product Name: Recombinant Human LARP4B Protein
Catalog #: PEH8828



Summary

Name	LAR4B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92615
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PADMNALALG PSEYDSLPEP SETGGNESQP DSQEDPREVL KKTLEFCLSR ENLASDMYLI SQMDSQYVP ITTVANLDHI KKLSTDVDLI VEVLRLPLV QVDEKGEKVR PNQNR CIVIL REISESTPVE EVEALFKGDN LPKFINCEFA YNDNWFITFE TEADAQQAYK YLREEVKTFQ GKPIKARIKA

image



Background

Alternative Names	La-related protein 4B (La ribonucleoprotein domain family member 4B) (La ribonucleoprotein domain family member 5) (La-related protein 5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LARP4B Protein
Catalog #: PEH8828

