


Product Name: Recombinant Human LURAP1 Protein
Catalog #: PEH8831



Summary

Name	LURA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96LR2
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLRDVEGKVG RKTPEGLLRG LRGECELGTS GALLPGASS TGHDLGDKIM ALKMELAYLR AIDVKILQQL VTLNEGIEAV RWLLEERGTL TSHCSSLTSS QYSLTGGSPG RSRRGSWDSL PDTSTTDRLD SVSIGSFLDT VAPSELDEQG PPGAPRSEMD WAKVIAGGER ARTEVDVAAT RLGSLRAVWK PPGERLQGGP

image



Background

Alternative Names Leucine rich adaptor protein 1 (Leucine repeat adapter protein 35A)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LURAP1 Protein
Catalog #: PEH8831

