


**Product Name: Recombinant Human LRRC1 Protein**  
**Catalog #: PEH8836**



## Summary

|                              |                                                                                                                                                                                                                 |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | LRRC1_HUMAN                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                             |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                           |
| <b>Accession #</b>           | Q9BTT6                                                                                                                                                                                                          |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                          |
| <b>Species</b>               | Human                                                                                                                                                                                                           |
| <b>Molecular Mass</b>        | 21 kDa                                                                                                                                                                                                          |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                       |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                              |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.  |
| <b>Sequence</b>              | CVRDNRRLTRI PAEVSQATEL HVLDVAGNRL LHLPLSLTAL KLKALWLSDN<br>QSQPLLTFQT DTDYTTGEKI LTCVLLPQLP SEPTCQENLP RCGALENLVN<br>DVSDEAWNER AVNRVSAIRF VEDEKDEEDN ETRTLLRRAT PHPGELKHKM<br>KTVENLRNDM NAAKGLDSNK NEVNHAIDRV |

## image



## Background

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Leucine-rich repeat-containing protein 1 (LANO adapter protein) (LAP and no PDZ protein) |
|--------------------------|------------------------------------------------------------------------------------------|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human LRRC1 Protein**  
**Catalog #: PEH8836**

---

