


Product Name: Recombinant Human L1TD1 Protein
Catalog #: PEH8837



Summary

Name	LITD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T7N2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFARLAKKKE NITYMKREQL TETDKDIAPV LDLKCKDVSA IMNKFVKVLM IQDLMFEEMR ETLKNDLKAV LGGKATIEPV KNSSENSSRT EFQQIINLAL QKTGMVVGKIE GENSKIGDDN ENLTFKLEVN ELSGKLDNTN EYNSNDGKKL PQGESRSYEV MGSMEETLCN IDDRDGNRNV HLEFTERESR KDGEDEFVKE MREERKFQKL KNKEEVLKAS

image



Background

Alternative Names	LINE-1 type transposase domain-containing protein 1 (ES cell-associated protein 11)
--------------------------	---

Note

Product Name: Recombinant Human L1TD1 Protein
Catalog #: PEH8837



For Research Use Only , Not for Diagnostic Use