


Summary

Name	LOXH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVV2
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GMARYHVTVC TGELEGAGTD ANVYLCLFGD VGDTGERLLY NCRNNTDLFE KGNADFTIE SVTMRNVRV RIRHDGKGS SGWYLDRVLV REEQPESDN VEFPCRWLD KDKDDGQLVR ELLPSDSSAT LKNFRYHISL KTGDVSGAST DSRVYIKLYG DKSDTIKQVL LVSDNNLKDY FERGRVDEFT LETLNIGNIN RLVIGHDSTG MHASWFLGSV QIRVPRQGKQ YTFPANRWLD KNQADGRLEV ELYPSEVVEI QKLVHYEVEI WTGDVGGAGT SARVYMQIYG EKGKTEVFLF SSRSKVFERA SKDTFQTDTF TIYIDLAL TKIRIRHDNT GNRAGWFLDR IDITDMNNEI TYYFPCQRWL AVEEDDGQLS RELLPVDESY VLPQSEEGRG GGDNNPLDNL

image



Background

Product Name: Recombinant Human LOXHD1 Protein
Catalog #: PEH8839



Alternative Names

Lipoxygenase homology domain-containing protein 1

Note

For Research Use Only , Not for Diagnostic Use