

Product Name: Recombinant Human PPFIA2 Protein
Catalog #: PEH8841



Summary

Name	LIPA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75334
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRESTSIHDM NDKLENELAN KEAILRQMEE KNRQLQERLE LAEQKLQQTM RKAETLPEVE AELAQRIAL TKAEERHGNI EERMRLHLEGQ LEEKNQELQR ARQREKMNEE HNKRLSDTVD RLLTESNERL QLHLKERMAA LEEKNVLIQE SETFRKNLEE SLHDKERLAE EIEKLRSELD QMKMRTGSLI EPTIPRTHLD TSAELRYSVG SLVDSQSDYR TTKVIRRR

image



Background

Alternative Names Liprin-alpha-2 (Protein tyrosine phosphatase receptor type f polypeptide-interacting protein alpha-2) (PTPRF-interacting protein alpha-2)

Note

Product Name: Recombinant Human PPFA2 Protein
Catalog #: PEH8841



For Research Use Only , Not for Diagnostic Use