


Product Name: Recombinant Human PPPIA3 Protein
Catalog #: PEH8842



Summary

Name	LIPA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75145
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDGRRGSALG PDEAGGELER LMVTMLTERE RLLETLREAQ DGLATAQLRL RELGHEKDSL QRQLSIALPQ EFAALTKELN LCREQLLERE EEIAELKAER NNTRLLEHL ECLVSRHERS LRMTVVKRQA QSPGGVSSEV EVLKALKSLF EHHKALDEKV RERLRMALER VAVLEEELEL SNQETLNLRE QLSRRRSGLE EPGKDG DGQT LANGLGPGGD SNRRTAEELE ALERQRAEVC QLRERLAVLC RQMSQLEEL GTAHRELKGA EEANSKLQRD LKEALAQRED MEERITTLEK RYLSAQREAT SLHDANDKLE NELASKESLY RQSEKSRQL AEWLDDAKQK LQQTLLQKAET LPEIEAQLAQ RVAALNKAEE RHGNFEERLR QLEAQLEEKN QELQRARQRE KMNDDHNKRL SETVDKLLSE SNERLQLHLK ERMGALEEKN SLSEEIANMK KLQDELLLNK EQLLAEMERM QMEIDQLRGR PPSSYSRSLP GSALELRYSQ APTLPSGAHL

image



Background

Product Name: Recombinant Human PPFIA3 Protein
Catalog #: PEH8842



Alternative Names

Liprin-alpha-3 (Protein tyrosine phosphatase receptor type f polypeptide-interacting protein alpha-3) (PTPRF-interacting protein alpha-3)

Note

For Research Use Only , Not for Diagnostic Use