


Product Name: Recombinant Human LIX1L Protein
Catalog #: PEH8844



Summary

Name	LIX1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVB5
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEALQEFWQM KQSRGADLKN GALVVYEMVP SNSPPYVCYV TLPGGSCFGS FQFCPTKAEA RRSAAKIALM NSVFNEHPSR RITDEFIEKS VSEALASFNG NREEADNPNT GIGAFRFMLE SNKGKSMLEF QELMTVFQLL HWNGSLKAMR ERQCSRQEV L AHYSHRALDD DIRHQMALDW VSREQSVPGA LSRELASTER ELDEARLAGK ELRFHKEKGD ILVLAAGQLG

image



Background

Alternative Names LIX1-like protein

Note

Product Name: Recombinant Human LIX1L Protein
Catalog #: PEH8844



For Research Use Only , Not for Diagnostic Use