


Product Name: Recombinant Human LONP2 Protein
Catalog #: PEH8846



Summary

Name	LONP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86WA8
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRKPKQDDDK RVIAIRPIRR ITHISGTLED EDEDEDNDDI VMLEKKIRTS SMPEQAHKVC VKEIKRLKKM

image



Background

Alternative Names	Lon protease homolog 2, peroxisomal (EC 3.4.21.53) (Lon protease-like protein 2) (Lon protease 2) (Peroxisomal Lon protease) (pLon)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use