


Product Name: Recombinant Human MPP2 Protein
Catalog #: PEH8849



Summary

Name	MPP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14168
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TFSNQPVPD AVRMVGIRKT AGEHLGVTFR VEGGELVIAR ILHGGMVAQQ GLLHVGDIIK EVNGQPVGSD PRALQELLRN ASGSVILKIL PSYQEPHLPR QVFKCHFDFY DPARDLIPC KEAGLRFNAG DLLQIVNQDD ANWWQACHVE GGSAGLIPSQ LLEEKRFV KRDLELTPNS GTLCGSLSGK KKKRMMYLTT KNAEFDRHEL LIYEEVARMP PFRKTLVLI GAQGVGRRSL KNKLIMWDPD RYGTTVPYTS RRPKDSEREG QGYFVSRGE MEADVGRAGRY LEHGEYEGNL YGTRIDSIRG VVAAGKVCVL DVNPQAVKVL RTAEFVPYVV FIEAPDFETL RAMNRAALES GISTKQLTEA DLRRTVEESS RIQRGYGHYF DLCLVNSNLE RTFRELQTAM

image



Background

Product Name: Recombinant Human MPP2 Protein
Catalog #: PEH8849



Alternative Names

MAGUK p55 subfamily member 2 (Discs large homolog 2) (Protein MPP2)

Note

For Research Use Only , Not for Diagnostic Use