



Product Name: Recombinant Human MFHAS1 Protein
Catalog #: PEH8850



Summary

Name	MFHA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4C4
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQAGLTEDRL QSALSYLHES GKLLYFEDSP ALKEHVFHNL TRLIDILNVF FQRDPSLLLH KLLLGTSGEG KAEGESSPPM ARSTPSQELL RATQLHQYVE GFLHGLLPA HVIRLLLKPH VQAQQLQLL LELLEKMGLC YCLNPKGKP LNGSTAWYKF PCYVQNEVPH AEAWINGTNL AGQSFVAEQL QIEYSFPFTF PLGLFARYSV QINSHVVHRS DGKFQIFAYR GKVPVVVSYP ARGVLPQPD LSIASHASLP NIWTAWQAIT PLVEELNVLL QEWPGHLHYTV HILCSKCLKR GSPNPHAFPG ELLSQRPPEG VAEIICPKNG

image



Background

Alternative Names	Malignant fibrous histiocytoma-amplified sequence 1 (Malignant fibrous histiocytoma-amplified sequence with leucine-rich tandem repeats 1)
--------------------------	--

Product Name: Recombinant Human MFHAS1 Protein
Catalog #: PEH8850



Note

For Research Use Only , Not for Diagnostic Use