


Product Name: Recombinant Human MAP6D1 Protein
Catalog #: PEH8851



Summary

Name	MA6D1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9H5
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQPPAGARDS GRDVPLTQYQ RDFGLWTPA GPKDPPPGRG PGAGGRRGKS SAQSSAPPAP GARGVYVLP I GDADAAA AVT TSYRQEFQAW TGVKPSRSTK TKPARVITTH

image



Background

Alternative Names MAP6 domain-containing protein 1 (21 kDa STOP-like protein) (SL21)

Note

For Research Use Only , Not for Diagnostic Use