



Product Name: Recombinant Human MREG Protein
Catalog #: PEH8854



Summary

Name	MREG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N565
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERALPEKEPL VSDNNPYSSF GATLVRDDEK NLWSPHVDVS HTEADDDRTL YNLIVIRNQQ AKDSEEWQKL NYDIHTLRQV RREVRNRWKC ILEDLGFQKE ADSLLSVTKL STISDSKNTR KAREMLLKLA EETNIFPTSW ELSERYLFVV DRLIALDAAE EFFKLARTY PCKPGVPCLA DGQKELHYLP

image



Background

Alternative Names Melanoregulin (Dilute suppressor protein homolog)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MREG Protein
Catalog #: PEH8854

