


Product Name: Recombinant Human MTX3 Protein
Catalog #: PEH8855



Summary

Name	MTX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5HYI7
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGGWGLPSVH SESLVVMAYA KFSGAPLKVN VIDNTWRGSR GDVPILTTED DMVSQPAKIL NFLRKQKYNA DYELSAKQGA DTLAYIALLE EKLLPAVLHT FWVESDNYFT VTKPWFASQI PFPLSLILPG RMSKGALNRI LLTRGQPPLY HLREVEAQIY RDAKECLNLL SNRLGTSQFF FGDTPSTLDA YVFGFLAPLY KVRFPKVQLQ EHLKQLSNLC RFCDDILSSY FRLSLGGISP AGQETVDANL

image



Background

Alternative Names Metaxin-3

Note

Product Name: Recombinant Human MTX3 Protein
Catalog #: PEH8855



For Research Use Only , Not for Diagnostic Use