


Product Name: Recombinant Human MYH13 Protein
Catalog #: PEH8868



Summary

Name	MYH13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKX3
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVKSELDKRV IEKDEEIEQL KRNSQRAAEA LQSVLDAEIR SRNDALRLKK KMEGDLNEME IQLGHSNRQM AETQKHLRTV QGQLKDSQLH LDDALRSNED LKEQLAIVER RNGLLLEELE EMKVALEQTE RTRRLSEQEL LDASDRVQLL HSQNTSLINT KKKLEADIAQ CQAEVENS IQ ESRNAEEKAK KAITDAAMMA EELKKEQDTS AHLERMKKNL EQTVKDLQHR LDEAEQLALK GGKKQIQKLE NRVRELENEL DVEQKRGAEA LKGAHKYERK VKEMTYQAE E DHKNILRLQD LVDKLQAKVK SYKRQAEAE EQANTQLSRC RRVQHELEEA AERADIAESQ VNKLRAKSRD

image



Background

Alternative Names	Myosin-13 (Myosin heavy chain 13) (Myosin heavy chain, skeletal muscle,
--------------------------	---

Product Name: Recombinant Human MYH13 Protein
Catalog #: PEH8868



extraocular) (MyHC-EO) (Myosin heavy chain, skeletal muscle, laryngeal)
(MyHC-IIL) (Superfast myosin)

Note

For Research Use Only , Not for Diagnostic Use