


Product Name: Recombinant Human MYH15 Protein
Catalog #: PEH8869



Summary

Name	MYH15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2K3
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDEEIENFRR KQQTIDSLQ SSLDSEAKSR IEVTRLKKKM EEDLNEMELQ LSCANRQVSE ATKSLGQLQI QIKDLQMQLD DSTQLNSDLK EQVAVAERRN SLLQSELEDL RSLQEQTERG RRLSEEELE ATERINLFYT QNTSLLSQKK KLEADVARMQ KEAEVVQEC QNAEEKAKKA AIEAANLSEE LKKKQDTIAH LERTRENMEQ TITDLQKRLA EAEQMALMGS RKQIQKLESR VRELEGELEG EIRRSAEAQR GARRLERCIK ELTYQAEEDK KNLSRMQTQM DKLQLKVQNY KQQVEVAETQ ANQYLSKYKK QQHELNEVKE

image



Background

Alternative Names Myosin-15 (Myosin heavy chain 15)

Product Name: Recombinant Human MYH15 Protein
Catalog #: PEH8869



Note

For Research Use Only , Not for Diagnostic Use