


Product Name: Recombinant Human NAA60 Protein
Catalog #: PEH8871



Summary

Name	NAA60_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H7X0
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDIDTVKHLG GDFWPIEYPD SWYRDITSNK KFFSLAATYR GAIVGMIVAE IKNRTKIHK E DGDILASNFS VDTQVAYILS LGVVKEFRKH GIGSLLES L KDHISTTAQD HCKAIYLHVL TTNNTAINFY ENRDFKQHHY LPYYYSIRGV LKDGFTYVLY

image



Background

Alternative Names	N-alpha-acetyltransferase 60 (hNaa60) (EC 2.3.1.259) (Histone acetyltransferase type B protein 4) (HAT4) (EC 2.3.1.48) (N-acetyltransferase 15) (N-alpha-acetyltransferase F) (NatF)
--------------------------	--

Note

Product Name: Recombinant Human NAA60 Protein
Catalog #: PEH8871



For Research Use Only , Not for Diagnostic Use