


Product Name: Recombinant Human NUDCD3 Protein
Catalog #: PEH8883



Summary

Name	NUDC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVD9
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	METGAAELYD QALLGILQHV GNVQDFLRVL FGFLYRKTDF YRLLRHPSDR MGFPPGAAQA LVLQVFKTFD HMARQDDEKR RQELEEKIRR KEEEEAKTVS AAAAEKEPVP VPVQEIEIDS TTELDGHQEV EKVQPPGPVK EMAHGSQEAE APGAVAGAAE VPREPPILPR IQEQFQKNPD SYNGAVRENY TWSQDYTDLE VRVPVPKHVV KKGQVSVALS SSSIRVAMLE ENGERVLMEG KLTHKINTES SLWSLEPGKC VLVNLSKVGE YWWNAILEGE EPIDIDKINK ERSMATVDEE EQAVLDRLTF DYHQKLQGKP QSHELKVHEM LKKGWDAEGS PFRGQRFDPA MFNISP GAVQ F

image



Background

Alternative Names NudC domain-containing protein 3

Product Name: Recombinant Human NUDCD3 Protein
Catalog #: PEH8883



Note

For Research Use Only , Not for Diagnostic Use