


Product Name: Recombinant Human OTUD1 Protein
Catalog #: PEH8884



Summary

Name	OTUD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VV17
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGPWGEEHLA ERGPRGWERG GDRCDAPGGD AARRPDPEAE APPAGSIEAA PSSAAEPVIV SRSDPRDEKL ALYLAEVEKQ DKYLRQRNKY RFHIIPDGNC LYRAVSKTVY GDQSLHRELK EQTVHYIADH LDHFSPLIEG DVGEFIIAAA QDGAWAGYPE LLAMGQMLNV NIHLTTGGRL ESPTVSTMIH YLGPEDSLRP SIWLSWLSNG HYDAVFDHSY PNPEYDNWCK QTQVQRKRDE ELAKSMAISL SKMYIEQNAC

image



Background

Alternative Names	OTU domain-containing protein 1 (EC 3.4.19.12) (DUBA-7)
--------------------------	---

Product Name: Recombinant Human OTUD1 Protein
Catalog #: PEH8884



Note

For Research Use Only , Not for Diagnostic Use