


Product Name: Recombinant Human OTUD3 Protein
Catalog #: PEH8885



Summary

Name	OTUD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T2D3
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PESGGGGGCE EEFVSFANQL QALGLKLREV PGDGNCLFRA LGDQLEGHSR NHLKHRQETV DYMIKQREDF EPFVEDDIPF EKHVASLAKP GTFAGNDAIV AFARNHQLNV VIHQLNAPLW QIRGTEKSSV RELHIAYRYG EHYDSVRRIN DNSEAPAHLQ TDFQMLHQDE SNKREIKTK GMDSEDDL RD EVEDAVQKVC NATGCSDFNL IVQNLEAENY NIESAIIAVL RMNQGKRNN A EENLEPSGRV LKQCGPLWEE GSGGARIFGN QGLNEGR TEN NKAQASPSEE NKANKNQLAK VTNKQRREQQ WMEKKKRQEE

image



Background

Alternative Names OTU domain-containing protein 3 (EC 3.4.19.12)

Product Name: Recombinant Human OTUD3 Protein
Catalog #: PEH8885



Note

For Research Use Only , Not for Diagnostic Use