


Product Name: Recombinant Human OSGIN1 Protein
Catalog #: PEH8886



Summary

Name	OSGI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJX0
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YTKPDAIHPH PLLQRKLTEA PGVSILDQDL DYLSEGLEGR SQSPVALLFD ALLRPDTDFG GNMKSVLTWK

image



Background

Alternative Names	Oxidative stress-induced growth inhibitor 1 (Bone marrow stromal cell-derived growth inhibitor) (BMSC-derived growth inhibitor) (Ovary, kidney and liver protein 38) (huOKL38) (Pregnancy-induced growth inhibitor OKL38)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use