


Product Name: Recombinant Human OSBPL10 Protein
Catalog #: PEH8887



Summary

Name	OSB10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXB5
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DNEDKEETEL GVMEDQRSII LHLISQLKLG MDLTKVVLPT FILEKRSLE MYADFMAHPD LLLAITAGAT PEERVICFVE YYLTAFFEGR KGALAKKPN PIIGETFHCS WEVPKDRVKP KRTASRSPAS CHEHPMADDP SKSYKLRFA EQVSHHPPIS CFYCEEEKR LCVNTHVWTK SKFMGMSVGV SMIGEGVLR LEHGEEYVFT LPSAYARSIL TIPWVELGGK VSINCAKTGY SATVIFHTKP FYGGKVHRVT AEVKHNPTNT IVCKAHGEWN GTLEFTYNGG ETKVIDTTTL PVYPKIRPL EKQGPMSRN LWREVTRYLR LGDIDAATEQ KRHLEEKQRV EERKRENLRT PWKPKYFIQE GDGWVYFNPL

image



Background

Alternative Names Oxysterol-binding protein-related protein 10 (ORP-10) (OSBP-related protein)

Product Name: Recombinant Human OSBPL10 Protein
Catalog #: PEH8887



10)

Note

For Research Use Only , Not for Diagnostic Use